

MEDICAL REVITALIZATION, INC Hill Rom P8000



Features

- New Mattress
- New Paint
- Pumps and Stickers

RefurbishmentProcess

- Disassemble Bed
- Steam clean frame and components
- Paint frame
- Install new brake mechanisms on both sides
- Sand and repaint shrouds with our special 2 part polyurethane coating

Specifications

- Maximum Weight Limit: 700lbs.
- Height: Low Position 20.7", High Position: 34.25"
- Overall Length: 83"
- Overall Width: 36", (Side Rails up: 34.25")
- Caster Size: 8"
- SideRail Length: 47"
- Maximum Head Elevation: 90°

- Rebuild pumps with new seals and adjust pumps speed and travel
- New stickers per customers choice of color
- Clean, adjust and lube brake/steer components
- New bumpers and plastic
- 1 year warranty