1-818-899-8090



Tuttnauer EZ10



Features:

- Fully automatic, one touch design
- Chamber door stays closed during the entire cycle
- Closed door active drying system to maintain sterility and efficient drying of packs and pouches
- 316L grade stainless steel with electro-polish finish
- Double safety locking device prevents door from opening while chamber is pressurized
- 4 trays
- 510k clearance for dental handpieces
- Drain valve located in the front, allowing for a quick and easy draining of water reservoir
- Customize and store program parameters

Specifications:

- Chamber Dimensions: 10" x 19"
- Chamber Volumes: 6gal
- Overall Dimensions: 21.50" x 20" x 14.40"
- Standard Cassette Capacity: 3 Full and 3 Half
- Standard Cycle Time Unwrapped: Cold 30min, Hot 14min
- Power Input: 120V, 50hz/60Hz