



MEDICAL REVITALIZATION, INC

Midmark 7300



Specifications:

- Height: 27.50" to 44"
- Height: 27.50" to 44"
- Trendelenburg and Reverse Trendelenburg: 0-28 degrees
- Left and Right Lateral Tilt: 0-18°
- Back Section: 28°-0-28° (normal), 80°-0-105° (reverse)
- Normal Position Coordinated: flex 220°, reflex 155°
- Reverse position Coordinated: flex 236°, reflex 124°
- Leg section: down 105°, up 80° for chair positioning
- Single Button Return to level
- Weight Capacity: 400lbs

Refurbishment Process

- Fully Disassembled Table
- Clean, Degrease, and Sandblast Table Parts
- Repaint All Painted Parts
- Inside Base Motor Painted Red To Show Work Detail
- Table top frame and outside base painted our specialty coating hammered bronze for a great look and durable finish.
- Hydraulic block and lines flushed
- Circuit board and wire harnesses inspected and tested.
- Radiolucent tops smooth and recoated.
- New stickers and emblems
- Caster re-greased and sealed
- New brake feet pads
- New 2" pad set
- 1 Year Warranty

Features

- Designed for maximum C-arm access and X-ray visibility and uses electro/hydraulic power to move the table quickly from one position to another.
- There is a almost infinite variety of positions with independent head, back, seat, and leg sections.
- The unique pad design prevents the patient from sinking inot the gap between the X-ray tops.
- If the patient's position is reversed, the micro-processor reverses command.
- The Midmark 7300 has many optional accessories available and also accepts other manufacturers accessories.